

# Choose Cost-effective Restek LC Kits

## For Your Preventative Maintenance



- Significant savings over instrument manufacturer prices.
- High-quality components in every kit.
- Wide range of options for HPLC systems and pumps.
- Full product listing on reverse side.



[www.restek.com/LC-Maintenance](http://www.restek.com/LC-Maintenance)



# LC Maintenance Kits


- Significant savings over instrument manufacturer prices.
- High quality components in every kit.
- Wide selection of PM Kits for your HPLC systems.

Description for Agilent LCs	Model #	Similar to Agilent part #	qty.	cat.#
Preventive Maintenance Kit Includes: rotor seal, needle seat, needle assembly, seat cap	1050	01078-68721	kit	25259
Autosampler Preventive Maintenance Kit Includes: rotor seal, needle assembly, needle seat	1100, 1200	G1313-68709	kit	25271
Pump Maintenance Kit Includes: PTFE frit, outlet cap, active inlet cartridge, gold disk seal, 2 piston seals, glass solvent filter	1050, 1100, 1200	G1311-68710	kit	25270
Description for Waters LCs	Model #	Similar to Waters part #	qty.	cat.#
Preventive Maintenance Kit Includes: sparge diffuser, filter insert, compression screws, SS ferrule, battery for 2690/717, 250µL WISP syringe, seal wash plunger seal kits (2), wash tube seal kits (4), 2690 seal pack rebuild kit (steel bodies not included), 2690 head plunger seal kits (2), solvent reservoir 20µm filters (4), Alliance check valve cartridges (2), Alliance plunger assemblies (2), 2690 face seals (4)	Alliance 2690, 2695	WAT270944	kit	25143
Preventive Maintenance Kit Includes: PerformancePLUS cartridges (2), sparge diffusers (4), Super Seals (2), solvent reservoir 20µm filters (2), sapphire plungers (2), reference valve rebuild kit, inlet manifold kit	600 Pump	WAT052675	kit	25144
Preventive Maintenance Kit Includes: 717 seal pack with needle, filter insert, 250µL WISP syringe	717 Autosampler	WAT052669	kit	25145
Preventive Maintenance Kit Includes: 616/326/625 plunger assemblies (2), pump seal kit, sparge diffusers (4), solvent reservoir 20µm filters (4), 616 cartridge assemblies (4)	616 Pump	WAT052672	kit	25146
Preventive Maintenance Kit Includes: PerformancePLUS check valves (2), sparge diffuser, solvent 10µm filter, sapphire plungers (2), plunger seals (2)	515 Pump	WAT052587	kit	26519
Preventive Maintenance Kit Includes: PerformancePlus cartridges (8), Clear-100 seals (4), sapphire plungers (4), solvent filters (2), reference valve rebuild kit	1525 Pump	201000114	kit	26430



To order, visit  
[www.restek.com/LC-Maintenance](http://www.restek.com/LC-Maintenance)

Restek also offers lamps, seals, syringes, and other LC accessories at [www.restek.com/LCacc](http://www.restek.com/LCacc)



#### PATENTS & TRADEMARKS

Restek patents and trademarks are the property of Restek Corporation. Other trademarks appearing in Restek literature or on its website are the property of their respective owners. The Restek registered trademarks used here are registered in the United States and may also be registered in other countries.



**Restek U.S.** • 110 Benner Circle • Bellefonte, PA 16823 • 814-353-1300 • 800-356-1688 • fax: 814-353-1309 • [www.restek.com](http://www.restek.com)  
**Restek France** • phone: +33 (0)1 60 78 32 10 • fax: +33 (0)1 60 78 70 90 • e-mail: [restek@restekfrance.fr](mailto:restek@restekfrance.fr)  
**Restek GmbH** • phone: +49 (0)6172 2797 0 • fax: +49 (0)6172 2797 77 • e-mail: [info@restekgmbh.de](mailto:info@restekgmbh.de)  
**Restek Ireland** • phone: +44 (0)2890 814576 • fax: +44 (0)2890 814576 • e-mail: [restekireland@aol.com](mailto:restekireland@aol.com)  
**Restek Japan** • phone: +81 (3)6459 0025 • fax: +81 (3)6459 0025 • e-mail: [restekjapan@restek.com](mailto:restekjapan@restek.com)  
**Thames Restek U.K. LTD** • phone: +44 (0)1494 563377 • fax: +44 (0)1494 564990 • e-mail: [sales@thamesrestek.co.uk](mailto:sales@thamesrestek.co.uk)

Lit. Cat.# GNTS1527-UNV

© 2012 Restek Corporation. All rights reserved.  
 Printed in the U.S.A.

